

**Product Name :**  
Naugra Lab Air Curtain

**Product Code :**  
NLAB-LABORATORYEM6603



**Description :**

Naugra Lab Air Curtain

**Technical Specification :**

Air Curtain

A must for every laboratory, clean rooms, industries etc.

All Air Curtains when installed on a door forms an invisible curtain/barrier of continuous air and thus prevents escape of conditioned/treated air and entry of contaminated hot, humid, dirty and unwanted polluted air inside the room.

Most effective insulation for energy conservation and preventing temperature loss in controlled atmospheric areas, thus maintaining very high levels of cleanliness by forming a very strong barrier to prevent movement of dust particles across it.

An air Curtain is a modern device which, when installed on a door, forms

an invisible curtain of continuous air and thus prevents escape of conditioned air and entry of outside hot, humid, dirty and unwanted polluted air.

The cabinets of Air Curtain are made of cold rolled mild steel sheets.

A must for every laboratory, clean rooms, industries etc.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003