

Product Name :
Nauga Lab Anderson Pipette

Product Code :
NLAB-LABORATORYEM65004



Description :

Nauga Lab Anderson Pipette

Technical Specification :

Anderson Pipette

Particle size determination apparatus, (sedimentation pipette, anderson) consisting of 500ml.

Graduated cylinder, graduation from 0 to 20cm and fitted with 10rni.

Pipette with three way stopcock, the tip of pipette at the level of 0 mark at the cylinder

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

