


Product Name :
Nauga Lab Antibiotic Zone Reader

Product Code :
NLAB-LABORATORYEM65003



Description :

Nauga Lab Antibiotic Zone Reader

Technical Specification :

Antibiotic Zone Reader

The Antibiotic zone reader provides a rapid accurate method for determining the strength of antibiotic materials by measuring the diameter of an inhibited zone in petri dish.

This zone reader measures the diameter of the inhibited zone to 0.02mm with the range of 0 to 35mm.

It works on 220V AC 50Hz

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com