

Product Name :
Naugra Lab Autoclave Laboratory Portable

Product Code :
NLAB-LABORATORYEM47003



Description :

Naugra Lab Autoclave Laboratory Portable

Technical Specification :

Autoclave Laboratory Portable

Made of heavy duty stainless steel

Double safety valve.

Fitted with 2.5" special Gauge

Steam release valve

Portable
Autoclaves dia. x Height in mm:

300 x 300mm (Seamless)

300x37mm (Joined)

300x50mm (Joined)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003