

Product Name :
Naugra Lab Ball Mill

Product Code :
NLAB-LABORATORYEM65007



Description :

Naugra Lab Ball Mill

Technical Specification :

Ball Mill

Consists of following features:-

Electrically operated, motor driven made of stainless steel.

Capacity : 1, 2 and 5 litre mill made of stainless steel

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

