

Product Name :
Naugra Lab B.O.D. Incubator

Product Code :
NLAB-LABORATORYEM33001



Description :

Naugra Lab B.O.D. Incubator

Technical Specification :

Double wall
construction.

Inner chamber made
richly anodized aluminium or highly polished stainless steel sheet.

Supplied with two/three
removable perforated shelves.

Exterior made out of thick
mild steel sheet duly finished in white stoving enamel/powder coated paint with
mat finished colour combination.

Outer double wall door
is provided with magnetic gasket and lock & key arrangement.

Inner door unbreakable acrylic transparent sheet for inspection of material in the chamber without disturbing the chamber temperature.

Optional:

Timer 0-24 hours day night cycle.

Co2 /air mixer nozzle is provided.

Internal illumination with 3 fluorescent tubes.

Microprocessor PID digital temp. Indicator-cum-controller.

Automatic Voltage stabilizer 3-5 KVA (recommended-must).

Capacity in Liter Size in mm WxDxH

205 x 405

452 x 610

675 x 830

900 x 865

620x580x900

720x680x900

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003