

**Product Name :**  
Naugra Lab B.O.D. Incubator

**Product Code :**  
NLAB-LABORATORYEM33001



**Description :**

Naugra Lab B.O.D. Incubator

**Technical Specification :**

Double wall construction.

Inner chamber made richly anodized aluminium or highly polished stainless steel sheet.

Supplied with two/three removable perforated shelves.

Exterior made out of thick mild steel sheet duly finished in white stoving enamel/powder coated paint with mat finished colour combination.

Outer double wall door is provided with magnetic gasket and lock & key arrangement.

Inner door unbreakable acrylic transparent sheet for inspection of material in the chamber without disturbing the chamber temperature.

**Optional:**

Timer 0-24 hours day night cycle.

Co2 /air mixer nozzle is provided.

Internal illumination with 3 fluorescent tubes.

Microprocessor PID digital temp. Indicator-cum-controller.

Automatic Voltage stabilizer 3-5 KVA (recommended-must).

**Capacity in Liter Size in mm WxDxH**

~~205 liter~~  
205 liter 305x405

~~452 liter~~  
452 liter 450x610

~~675 liter~~  
675 liter 675x830

~~960 liter~~  
960 liter 960x865

620x650x900

700x650x900

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003