

Product Name :
Naugra Lab Boss Heads Stands & Fixing Devices

Product Code :
NLAB-LABORATORYEM77016

Description :

Naugra Lab Boss Heads Stands & Fixing Devices

Technical Specification :

Boss Heads, Stands & Fixing Devices

All stands have perfectly machined and stable bases and machined perfectly vertical uprights with standard 9 mm diameter.

Boss Heads (x-blocks) are perfectly machined out of solid Brass rods.

Complete with 38 cm high upright.

Complete with 38 cm high upright with quick tangential adjustment & clamp nut.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003