

**Product Name :**  
Nauga Lab Bottle Cap Sealing Machine

**Product Code :**  
NLAB-LABORATORYEM65005



**Description :**

Nauga Lab Bottle Cap Sealing Machine

**Technical Specification :**

Bottle Cap Sealing Machine:-

Portable, hand operated equipment with two extra cap holders of different sizes.

Totally casted machine hand operated.

Hand operated.

Electrical operated (mounted on heavy stand)

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments