

Product Name :
Nauga Lab Bottle Cap Sealing Machine

Product Code :
NLAB-LABORATORYEM65005



Description :

Nauga Lab Bottle Cap Sealing Machine

Technical Specification :

Bottle Cap Sealing Machine:-

Portable, hand operated equipment with two extra cap holders of different sizes.

Totally casted machine hand operated.

Hand operated.

Electrical operated (mounted on heavy stand)

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments