

Product Name :
Nauga Lab Bottle Filling Machine

Product Code :
NLAB-LABORATORYEM76015



Description :

Nauga Lab Bottle Filling Machine

Technical Specification :

Bottle Filling Machine

Semi Automatic electric operated having filling cap. 50 to 300ml.

Complete stainless steel body.

With single syringe

With double syringe

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003