

Product Name :
Nauga Lab Bottle / Liquid Filling Machine

Product Code :
NLAB-LABORATORYEM65009



Description :

Nauga Lab Bottle / Liquid Filling Machine

Technical Specification :

Bottle / Liquid Filling Machine:-

A sturdy and heavy unit made of CRC sheet mounted on a heavy angle iron duly powder coated.

Capacity of output of water is adjustable through cam system knob provided from 50 ml. to 300 ml.

It can be operated in two ways:-

Hand operated.

Motorized (Semi Automatic) operated

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003