

**Product Name :**  
Nauga Lab Bottle / Liquid Filling Machine

**Product Code :**  
NLAB-LABORATORYEM65009



**Description :**

Nauga Lab Bottle / Liquid Filling Machine

**Technical Specification :**

Bottle / Liquid Filling Machine:-

A sturdy and heavy unit made of CRC sheet mounted on a heavy angle iron duly powder coated.

Capacity of output of water is adjustable through cam system knob provided from 50 ml. to 300 ml.

It can be operated in two ways:-

Hand operated.

Motorized (Semi Automatic) operated

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003