

Product Name :
Naugra Lab Bottle Washing Machine

Product Code :
NLAB-LABORATORYEM65008



Description :

Naugra Lab Bottle Washing Machine

Technical Specification :

Bottle Washing Machine:-

Washing tank made of stainless steel mounted on sturdy square ms pipe frame duly powder coated.

Motorized driven unit embedded in a metallic box rested on the main unit with a drill check to accommodate various size of washing brushes.

Workable on 220V-50Hz AC mains

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

