

Product Name : Naugra Lab Cannula Borosil Glass	Product Code : NLAB-LABORATORYEM65013
-----------------------------------------------------------	-------------------------------------------------

--

Description :

Naugra Lab Cannula Borosil Glass

Technical Specification :

Cannula Borosil Glass

Made from Borosilicate Glass:-

Arterial

Venous

Syme's Perfusion

T-Shaped

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003