


Product Name :
Naugra Lab Capsule Filling Machine

Product Code :
NLAB-LABORATORYEM65012



Description :

Naugra Lab Capsule Filling Machine

Technical Specification :

Capsule Filling Machine

A compact hand operated machine suitable to fill two-part hard gelatine capsules no. 0, 1 & 2 by interchangeable plates.

Main body made of heavy cast iron and all contact parts made of brass richly plated or stainless steel.

Handle locking arrangement makes the operator to lift the loading plate without fear of damaging the capsules.

Pressing is achieved by cam system which also saves time to press the plunger as compared to screw system.

Plate of stainless steel for any one no. Is supplied with machine.

For filling 100 to 200 capsules.

Optional:- Plates of other nos. supplied, on demand

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003