

**Product Name :**  
Nauga Lab Centrifuge Angle Rotor Head Doctor Pattern

**Product Code :**  
NLAB-LABORATORYEM47004



**Description :**

Nauga Lab Centrifuge Angle Rotor Head Doctor Pattern

**Technical Specification :**

Centrifuge Angle Rotor Head Doctor Pattern

4000 to 5000 r.p.m. with speed regulator and plug & cord.

Capacity:

4x15ml tubes

6x15ml tubes
8x15ml tubes
12x15ml tubes

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003