

Product Name :
Nauga Lab Centrifuge Angle Rotor Head Doctor Pattern

Product Code :
NLAB-LABORATORYEM47004



Description :

Nauga Lab Centrifuge Angle Rotor Head Doctor Pattern

Technical Specification :

Centrifuge Angle Rotor Head Doctor Pattern

4000 to 5000 r.p.m. with speed regulator and plug & cord.

Capacity:

4x15ml tubes

6x15ml tubes
8x15ml tubes
12x15ml tubes

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003