

Product Name :
Nauga Lab Centrifuge Rectangular

Product Code :
NLAB-LABORATORYEM48002



Description :

Nauga Lab Centrifuge Rectangular

Technical Specification :

Centrifuge Rectangular

4000 to 5000 R.P.M.
fitted with electronic solid state speed regulator.

Top cover made of
transparent acrylic sheet, cord & plug. Optional:

Automatic Timer (Manual)

Automatic Timer
(Digital)

Speed Meter (Manual)

Speed Meter (Digital)

Capacity
4x15ml tubes.
6x15ml tubes.
8x15ml tubes.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003