

Product Name :
Naugra Lab Centrifuge

Product Code :
NLAB-LABORATORYEM48001



Description :

Naugra Lab Centrifuge

Technical Specification :

Electrically operated, dome type, sturdy, stable heavy body with resiliently mounted motor working on 220V AC 50 Hz single phase for vibration free performance with built in 5-speed regulator and a non-removable swing out rotor.

Complete with tube carriers.

Maximum speed 3500 RPM.

Neatly finished in staving enamel paint.

Sizes:-

4x15ml tubes.
6x15ml tubes.
8x15ml tubes.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003