

Product Name :
Naugra Lab CNC Lathe System

Product Code :
NLAB-EDUCNCLAB110036



Description :

Naugra Lab CNC Lathe System

Technical Specification :

Distance between centers: 600mm

Spindle bore: ?26mm

Spindle speeds: 30-4000rpm

Pitch-metric : 0.2 – 3.5mm

Pitch-inch : 8-56TPI

X travel: 130mm

X repeatability positioning accuracy: 0.012

Z repeatability positioning accuracy: 0.013

Z travel: 465 mm

X rapid speed : 3000mm/min

Turning ? via lathe bed: ?260mm

Z rapid speed : 2500mm/min

X positioning accuracy: 0.03

Z positioning accuracy: 0.032

Motor power: 2.2KW 230V~50Hz

Packing size : 2220 x 820 x 1720mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003