

Product Name :
Nauga Lab CNC Milling Machine with Cabinet and PC

Product Code :
NLAB-EDUCNCLAB110037



Description :

Nauga Lab CNC Milling Machine with Cabinet and PC

Technical Specification :

Stepping Resolution 0.003mm

Mechanical repeatability 0.01mm

Spindle Speed 1100 - 11000 RPM

Travel: X = 300mm, Y = 146mm, Z = 153mm

Adjustable tapered brass gib on X-axis.

Ground steel sliding ways with adjustable brass gib.

Hard anodized aluminum table.

Rapid Travel of 30 in/min.

Universal Standard G-Code import.

1/2" (12.7mm) 20 pitch Cr-Mo leadscrews with adjustable split bronze nuts on all axes.

It has coolant system with the removable collector, mechanical storage system on the bottom of a Computer system, Controller and for Electrical wiring.

The system supplied with PC with Intel Pentium-IV CPU @ 2.4 GHz, 512MB RAM, 52X DVD Writer, Optical Mouse, 15" TFT Monitor & Keyboard

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003