

Product Name :
Nauga Lab CNC Milling Machine

Product Code :
NLAB-EDUCNCLAB110035



Description :

Nauga Lab CNC Milling Machine

Technical Specification :

We are committed to offering the most advanced CNC milling with comprehensive services and technical support.

Our technology is for mould/dies industry.

Our CNC Milling machine has high performance and high-speed cutting.

Features:-

X, Y, Z-axis travel: 900/600/500mm

Built-in spindle: 24000rpm

Max. Table load: Unrestricted

Cutting Feedrate: 1 ~ 20,000mm/min

Automatic Tool Changer: 20 Tools

Machine Type:-

CNC Milling Machine.

CNC Cutting Machine.

CNC Drilling Machine

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003