

Product Name :
Naugra Lab Conical Percolator

Product Code :
NLAB-LABORATORYEM65010



Description :

Naugra Lab Conical Percolator

Technical Specification :

Conical Percolator:-

Stainless steel jar on stand.

Capacity:- 500 ml, 1 lit. and 2 lit.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003