

**Product Name :**  
Naugra Lab Cyclo Mixer

**Product Code :**  
NLAB-LABORATORYEM76008



**Description :**

Naugra Lab Cyclo Mixer

**Technical Specification :**

Cyclo Mixer

These are bench top optical shakers incorporating platform for handling large tube racks, small flasks or reaction trays.

The shaker platform rotates in a horizontal plane and allows a wide variety of rotating and mixing applications.

These shakers have a built in 0-15 minutes digital timer.

Have a fix speed of 180 RPM and it is useful for sergological, VDRL tests and other tuned clinical diagnostics procedures requiring gentle agitation.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003