

Product Name :
Naugra Lab Cyclo Mixer

Product Code :
NLAB-LABORATORYEM76008



Description :

Naugra Lab Cyclo Mixer

Technical Specification :

Cyclo Mixer

These are bench top optical shakers incorporating platform for handling large tube racks, small flasks or reaction trays.

The shaker platform rotates in a horizontal plane and allows a wide variety of rotating and mixing applications.

These shakers have a built in 0-15 minutes digital timer.

Have a fix speed of 180 RPM and it is useful for sergological, VDRL tests and other tuned clinical diagnostics procedures requiring gentle agitation.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003