

**Product Name :**  
Naugra Lab Differential Blood Cell Counter

**Product Code :**  
NLAB-LABORATORYEM14004



**Description :**

Naugra Lab Differential Blood Cell Counter

**Technical Specification :**

Facility to count "8" Types of CELLS

1.16 Digits Alphanumeric LCD Display

Automatic percentage (%) calculation

Individual cells counting facility

In-built Stop Watch and Timer function

Audible beep facility

Differential Blood Cell Counter use microprocessor controlled system which offers accurate and quick differential counts.

This cell counter is featured with reliability, accuracy, stability

and is a sleek device.

These cell counters are available with 11 operating keys that are suitable for clinical labs for showing the primary bacterial or virus infections.

Specifications:-

**Cell Counting**  
**Key**

**Total 11 Cell**  
**keys and 3 Function keys)**

**Body**

**Display LCD**  
**Display**

**Power Input**  
**230V ± 10% AC,**  
**50 Hz**

Weight  
(Approx.)

Dimensions  
150 mm (LxBxH)

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003