

Product Name :
Nauga Lab Differential Blood Cell Counter

Product Code :
NLAB-LABORATORYEM14004



Description :

Nauga Lab Differential Blood Cell Counter

Technical Specification :

Facility to count "8" Types of CELLS

1.16 Digits Alphanumeric LCD Display

Automatic percentage (%) calculation

Individual cells counting facility

In-built Stop Watch and Timer function

Audible beep facility

Differential Blood Cell Counter use microprocessor controlled system which offers accurate and quick differential counts.

This cell counter is featured with reliability, accuracy, stability

and is a sleek device.

These cell counters are available with 11 operating keys that are suitable for clinical labs for showing the primary bacterial or virus infections.

Specifications:-

Cell Counting
Key

Total 11 Cell
keys and 3 Function keys)

Body

2 Line LCD
Display

Power Input
AC,
50 Hz

~~Weight~~
(Approx.)

~~Dimensions~~
150 mm (LxBxH)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003