

Product Name :
Nauga Lab Digital Electrophoresis Power Supply

Product Code :
NLAB-LABORATORYEM75002



Description :

Nauga Lab Digital Electrophoresis Power Supply

Technical Specification :

Digital Electrophoresis Power Supply

Delivers a current of any value between 1- 100 mA (500 Volts Max.), in a constant current mode.

Delivers a voltage of 0-500 Volts (100mA Max.), in constant voltage mode.

The range control potentiometer makes the output adjustable over the entire range of the unit.

Fitted with separate meters for voltage and current.

It has line and load regulation +0.5 % and ripples less than 10 mV.

All supplies are short circuit and overload protected.

Range 0-500 volts and 0-100 mA (With Two Digital Meters.)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003