

Product Name :
Naugra Lab Digital Fluorometer

Product Code :
NLAB-LABORATORYEM57001



Description :

Naugra Lab Digital Fluorometer

Technical Specification :

Digital Fluorometer

Features:-

3 Digit LED display

Covers entire visible range

Determines the concentration of Fluorophors & Quinine

Specifications:-

Fluorescence obtained with 1 ppm quinine sulphate in 0.1 N sulphuric acid (primary wavelength 360 nm)

Digital bright red seven segment LED display

2V, 50W Tungsten Halogen Lamp

Adjustable Ranges

Primary 5840 and Corning 5113 mounted on holders

Secondary 4378, Corning 3486, Corning 3385 mounted on holders

Highly sensitive wide range photodiode/cell

Sample Test Tube (H)

Minimum Sample
Power $\pm 10\%$, 50 Hz. AC
Dimensions 200 x 335 x 160 mm (Approx.)
Weight (Approx
Standard Accessories Secondary filters: A set of 5, Test Tube: A set of 5

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003