

**Product Name :**  
Nauga Lab Digital HB Reader LCD Alpha Numeric

**Product Code :**  
NLAB-LABORATORYEM76004



**Description :**

Nauga Lab Digital HB Reader LCD Alpha Numeric

**Technical Specification :**

Digital HB Reader LCD Alpha Numeric

Programmable with Memory

A very light weight, accurate & reliable Haemoglobin reader based on "Cyanmethemoglobin method".

It gives direct reading of unknown solution in concentration (gm/dl) without any calculations.

A microcontroller has been incorporated in the system for accuracy & reproducibility of the findings.

It operates in "Measurement & Calibration" mode.

The 15 gm/dl. Concentration of standard solution has been fed as a default value, however you can change this value by going into "Calibration" mode.

#### Technical Specifications:-

Wave Length Range:-	
---------------------	--

Minimum Volume:-	
------------------	--

Display	Alpha numeric
---------	---------------

Key Pad	Touch
---------	-------

--	--

<b>Light Source:-</b> LED	
<b>Detector:-</b> photo diode	
<b>Me Test:-</b> Storage Facility	
<b>Avail Zr:-</b>	
<b>Input supply:-</b> 10%, 50Hz, 1 AE	
<b>Battery:-</b>	
<b>Dimensions:-</b> 150 X 150 X 100 mm	

Test Tubes	10 X-10 mm (Plastic)
------------	----------------------

Weight- Net	
-------------	--

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003