

Product Name :
Naugra Lab Digital HB Reader LCD Alpha Numeric

Product Code :
NLAB-LABORATORYEM76004



Description :

Naugra Lab Digital HB Reader LCD Alpha Numeric

Technical Specification :

Digital HB Reader LCD Alpha Numeric

Programmable with Memory

A very light weight, accurate & reliable Haemoglobin reader based on "Cyanmethemoglobin method".

It gives direct reading of unknown solution in concentration (gm/dl) without any calculations.

A microcontroller has been incorporated in the system for accuracy & reproducibility of the findings.

It operates in "Measurement & Calibration" mode.

The 15 gm/dl. Concentration of standard solution has been fed as a default value, however you can change this value by going into "Calibration" mode.

Technical Specifications:-

Wavelength Range:-	
---------------------------	--

Minimum Volume:-	
-------------------------	--

Display:-	CD Alpha numeric
------------------	------------------

Key Pad:-	Touch
------------------	-------

--	--

Light Source -LED	
--------------------------	--

Detector -photo diode	
------------------------------	--

Meas. Facility -Storage Facility	
---	--

Auxiliary -	
--------------------	--

Main supply -10%, 50Hz, 1 AE	
-------------------------------------	--

Battery -	
------------------	--

Dimensions -100 mm	
---------------------------	--

Serial No	
Particulars	
Test Tubes-10 mm (Plastic)	

Weight	
- Net	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003