

Product Name :
Naugra Lab Digital Moisture Meter

Product Code :
NLAB-LABORATORYEM51002

Description :

Naugra Lab Digital Moisture Meter

Technical Specification :

Digital Moisture Meter

We are engaged in offering a wide range of Digital Moisture Meters.

Several products are offered under this product range that includes Universal Moisture Testing Machine, Digital Grain Testing Moisture Machine and Handy Digital Moisture Testing Meter.

Measuring Range:-

5% to 40%

Temperature Range:-

5 -50 Degree Celsius

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com