

**Product Name :**  
Naugra Lab Digital Moisture Testing Machine

**Product Code :**  
NLAB-LABORATORYEM51001



**Description :**

Naugra Lab Digital Moisture Testing Machine

**Technical Specification :**

**Digital Moisture Testing Machine**

Micro controller based  
Universal type, with rugged mechanical construction, accurate moisture content is direct readable on 3 digit L.E.D. display.

Temperature and  
conductive result are automatically compensated for final moisture percentage reading.

No conversion table is  
need.

The whole instruments is  
protected in wooden box.

**Features:-**

Easy operation, press and take reading

Direct Reading on digital display

No Adjustment

Results Consistent

**Specification:-**

Range :	
Accuracy :	
Dimension x 19 cm	
Weight : without box.	

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003