

Product Name :
Nauga Lab Digital Photo Fluorometer

Product Code :
NLAB-LABORATORYEM28001



Description :

Nauga Lab Digital Photo Fluorometer

Technical Specification :

Digital Photo Fluorometer

Specifications:-

Fluoro Sensitivity: Full-scale deflection obtained with 1 ppm quinine sulfate in 0.1 N sulphuric acid (primary wavelength 360 nm)

Display: 3 digits bright red seven segments LED display

Excitation Source: 12V, 50W Tungsten Halogen Lamp

Sensitivity Range: Adjustable in 4 ranges

Detector: Highly sensitive wide range photodiode/cell

Sample Test Tube Size: 15 mm (D) x 80 mm (H)

Minimum Sample: 4 ml

Power: 230 V \pm 10%, 50 Hz. AC

Dimension: 260 x 335 x 160 mm (Approx.)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003