

Product Name :
Naugra Lab Digital Photocolarimeter

Product Code :
NLAB-LABORATORYEM76002



Description :

Naugra Lab Digital Photocolarimeter

Technical Specification :

Digital Photocolarimeter

Specifications:-

Filters:- 5 filters 400 (V), 490 (B), 520 (G), 590 (Y), 650 (R) nm.

8 filters 40 (400), 49 (490), 52 (520), 54 (540),

Minimum Volume:- 1.0 ml.

Light Source:- 6 V, 10W, Tungsten Lamp.

Meter Scale:- A 135 mm. meter for direct measurement of optical density or % Transmission.

Photometric Accuracy:- + 1.0% of full scale.

Sensor:- Highly sensitive Photocell.

Power Consumption 20 Watts.

Weight Approx.:- 5.5 Kg

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003