

Product Name :
Nauga Lab Digital Photoelectric Colorimeter

Product Code :
PTH0003



Description :

Nauga Lab Digital Photoelectric Colorimeter

Technical Specification :

Digital Photoelectric Colorimeter

Specifications:-

Filters:- 8 filters 40 (400), 49 (490), 52 (520), 54 (540),

59 (590), 62 (620), 65 (650), 70 (700) nm.

Minimum Volume:- 1.0 ml.

Light Source:- 6 V, 10W, Tungsten Lamp.

Digital Display:- A 135 mm. meter for direct measurement of optical density or % Transmission.

Sensor:- Highly sensitive Photocell.

Power Consumption:- 230 V, 50 Hz. AC

Weight Approx.:- 5.5 Kg.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-0171-2643080, 2601773