

**Product Name :**  
Nauga Lab Digital Tele Thermometer

**Product Code :**  
NLAB-LABORATORYEM62001



**Description :**

Nauga Lab Digital Tele Thermometer

**Technical Specification :**

Digital Tele Thermometer

**Features:-**

Extremely useful for Pharmaceutical Labs.

Highly Accurate and Easy to Operate.

3½ Digit 7-Segment Display.

Six Temperature Probes

Digital Tele-Thermometer is an ideal, low cost, accurate and easy to operate instrument useful for pyrogen testing in pharmaceutical industries.

The instrument has 6 channels so that temperature of 6 different rabbits can be monitored simultaneously.

Special temperature probes are provided which can be easily inserted in the anus of the rabbits.

Any of the 6 channels can be selected by a rotary switch provided on the front panel.

**Specifications:-**

Range	0°C
Resolution	0.1°C
Accuracy	± 0.5°C, ± 1 Digit
RTD	(PT-100) or IC Sensor
Channels	6 Channels, switch selectable
Repeatability	± 0.1°C

Display, 7-segment bright red LED
Power ± 10% AC, 50 Hz.
Dimensions 75 mm (approx.)
Weight (approx.)
Accessories (Set of 6), Operation Manual, Dust Cover

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003