

Product Name :
Nauga Lab Disintegration Test Apparatus

Product Code :
NLAB-LABORATORYEM59002



Description :

Nauga Lab Disintegration Test Apparatus

Technical Specification :

Disintegration Test Apparatus

Features:-

Digital timer with 1 second resolution

Accurate time out provided by microprocessor based time totalizer

Baskets move Up & Down as per IP/USP test specifications

Temperature Controlled by sensitive thermostats

Digital Tablet Disintegration Apparatus is designed for accurate estimation of disintegration time of tablets as per IP/USP standards.

The instrument is designed to test to batches of six tablets, simultaneously.

The unit consists of two beakers having capacity of 1000 cc each.

In these beakers, suitable disintegration media can be taken.

These can be temperature controlled using a in-built heater and a highly thermostat fitted inside the instrument.

Microprocessor based time to talizer is provided for accurate time readout.

Specifications:-

Digital Tablet Disintegration Apparatus
30 and 20 rpm Speed
2% Accuracy
Resolution
Operating Temp. Range

Power 10%, 50Hz.AC
Dimensions (H) X 300 mm (W) X 320 mm (L) (Approx.)
Weight (Approx.)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003