

**Product Name :**  
Nauga Lab Dispensing Sampling Booth

**Product Code :**  
NLAB-LABORATORYEM44002



**Description :**

Nauga Lab Dispensing Sampling Booth

**Technical Specification :**

Dispensing Sampling Booth

Dispensing booth is a kind of partial purifying equipment for filling, refilling, weighing and sampling of raw material and compounds.

It is provided with HEPA filter, which prevents the airborne dusts by down draught technique.

The unit eliminates powder contamination to protect the operator and the surrounding environment.

Technical Specifications

Noise Level : Less than 65 dB

Vibration Level : Minimum

Power Supply : 230V AC,50HZ

MOC : GI Powder Coated / SS-304 / SS-316

Inside : Side panel with all three side with coving

Three Stage filtration

Pre-filters (90% down to 10 micron)

EU7 - Intermediate Filter (95% down to 3 micron)

EU14 - HEPA Filter - Supply / Exhaust (99.999% down to 0.3 micron)

Statically & Dynamically balanced Motor-blower with isolator to reduce noise level

Fluorescent Lights with Milky White Diffuser

Differential Pressure Gauges: 02 Nos. (Intermediate & HEPA Filter)

ON/OFF Switches (BL / FL / UV)

Auxiliary Powder Socket : 5/15 Amp for external equipment

DOP Port



|                         |  |
|-------------------------|--|
|                         |  |
| 11855<br>11855<br>11855 |  |
| 13890<br>13890<br>13890 |  |
| 13520<br>13520<br>13520 |  |
| 20880<br>20880<br>20880 |  |

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003