

Product Name :
Naugra Lab Distillation Apparatus All Glass

Product Code :
NLAB-LABORATORYEM71005



Description :

Naugra Lab Distillation Apparatus All Glass

Technical Specification :

Distillation Apparatus All Glass

The assembly consists of a corning glass flask with embedded heating elements fused in spiral type glass coil internally at bottom.

A double wall coiled condenser has standard ground glass joints.

Complete with glass fitting, metallic stand, rings clamps and electrical fittings.

Distillation Apparatus unit is as per NPL patent design.

Supplied with or without LOW WATER LEVEL CUT-OFF DEVICE (as ordered) for general safety of the heating element.

To work on 220/230 volts A.C. Supply.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003