

**Product Name :**  
Naugra Lab Double Cone Blender

**Product Code :**  
NLAB-LABORATORYEM65014



**Description :**

Naugra Lab Double Cone Blender

**Technical Specification :**

Double Cone Blender

It consists of two stainless steel cones welded with a cylindrical shell and having charging and discharging at opposite ends mounted on a channel stand.

Drive is by means of an electric motor & reduction gear box through a centre shaft.

The cone is mounted on plummer blocks for easy movement there by ensuring thorough mixing.

To facilitate easy movement of the cone while charging and discharging of the large units.

A special inching arrangement is provided with the unit motor used are single phase.

Capacity:-(kg)

5 kg.

10 kg

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003