

**Product Name :**  
Naugra Lab Dry Bath Incubator

**Product Code :**  
NLAB-LABORATORYEM72001



**Description :**

Naugra Lab Dry Bath Incubator

**Technical Specification :**

Dry Bath Incubator

Features :-

Designed for  
multi-function dry bath heating.

Ideal for Enzyme  
Reactions, Inactivation of Sera, Incubation and other laboratory procedures.

Precise  
temperature control from ambient 5OC to 100OC with 1OC accuracy.

User's self

compensation functions for temperature to counter for difference between block temperature and sample temperature.

Fitted with  
Microprocessor Based Pid Digital Temperature Indicator-cum-Controller for auto control with respect to time and temperature and is also fitted with In-Built Buzzer Timer For Auto Control To Indicate Cycle End.

Microprocessor  
Based Digital Temperature Indicator cum-Controller having individual display for set value and process value.

Heating cavity  
of 49mm carries one heating blocks of size 96mm (L) x 96mm (B) x 49 mm(H) each.

Supplied with  
any one set of heating block of various configuration available.

Recessed well/  
cavity accommodates one interchangeable aluminium insert block.

Block extraction  
tool facilitates lifting of hot blocks for easy interchangeability.

Dry Heating  
Block Incubators are powerful, compact and accurate everything you need for the precision heating of test tubes, cuvettes, micro centrifuge tubes and other small containers.

High performance  
and user friendly controls combine to produce a specification which offers exceptional efficiency and reliability.

The benefit to  
you-convenient, no fuss heating!

Excellent  
temperature control.

Uniform and  
stable heating

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003