

Product Name :
Naugra Lab Dry Bath Incubator

Product Code :
NLAB-LABORATORYEM72001



Description :

Naugra Lab Dry Bath Incubator

Technical Specification :

Dry Bath Incubator

Features :-

Designed for
multi-function dry bath heating.

Ideal for Enzyme
Reactions, Inactivation of Sera, Incubation and other laboratory procedures.

Precise
temperature control from ambient 5OC to 100OC with 1OC accuracy.

User's self

compensation functions for temperature to counter for difference between block temperature and sample temperature.

Fitted with
Microprocessor Based Pid Digital Temperature Indicator-cum-Controller for auto control with respect to time and temperature and is also fitted with In-Built Buzzer Timer For Auto Control To Indicate Cycle End.

Microprocessor
Based Digital Temperature Indicator cum-Controller having individual display for set value and process value.

Heating cavity
of 49mm carries one heating blocks of size 96mm (L) x 96mm (B) x 49 mm(H) each.

Supplied with
any one set of heating block of various configuration available.

Recessed well/
cavity accommodates one interchangeable aluminium insert block.

Block extraction
tool facilitates lifting of hot blocks for easy interchangeability.

Dry Heating
Block Incubators are powerful, compact and accurate everything you need for the precision heating of test tubes, cuvettes, micro centrifuge tubes and other small containers.

High performance
and user friendly controls combine to produce a specification which offers exceptional efficiency and reliability.

The benefit to
you-convenient, no fuss heating!

Excellent
temperature control.

Uniform and
stable heating

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003