

Product Name :
Nauga Lab Drying Cabinet/Oven

Product Code :
NLAB-LABORATORYEM40004



Description :

Nauga Lab Drying Cabinet/Oven

Technical Specification :

Drying Cabinet/Oven

This unit has double walled chamber with interior made of Stainless Steel and Exterior is either made of Mild Steel finished.

Standard Fitting:-

S.S Enclosure, Motor, Controller PID, Power Cord, Tube Heater, Rubber Gasket, S.S Tray, Door Lock.

Features:-

Designed, fabricated to suit various applications in Laboratory/Hospitals.

Double walled chamber with interior made of Stainless Steel and Exterior is either made of Stainless Steel/Mild Steel powder coated.

Gap of inner & outer walled filled with sufficiently thick Glass wool insulation.

Long lasting S.S. tubular heater are mounted at bottom to give uniform heating.

Motorized blower for maintaining uniform temp. with in chamber, blowing in Air from bottom.

Microprocessor based Controller cum Indicator with pt. 100 sensor.

Optional – P.I.D. Controller with Printer Interface P.I.D. Controller with PC Interface.

Sizes:-

| Capacity Size (Inch) | |
|-------------------------|--|
| 300 Ltr 300 x 300 | |
| 600 Ltr 400 x 400 | |
| 1500 Ltr 500 x 500 | |

| | |
|-----------------|--|
| 220 kVA approx | |
| 625 kVA approx | |
| 2000 kVA approx | |

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003