

Product Name :
Nauga Lab Drying Oven Industrial

Product Code :
NLAB-LABORATORYEM36001



Description :

Nauga Lab Drying Oven Industrial

Technical Specification :

Double walled and double door.

Temperature controlled by thermostat with dial type thermometer.

Forced air blower.

Timer 0-24 hours

Digital temp. Indicator

Size in mm Size in inches

900x600x600 (36"x24"x24")

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003