


Product Name :
Naugra Lab Egg Incubator

Product Code :
NLAB-LABORATORYEM34001



Description :

Naugra Lab Egg Incubator

Technical Specification :

Outer wall is made of
M.S. Sheet duly finished in white staving enamel/powder coated paint with mat
finished colour combination and inside stainless steel.

The gap between inner
and outer wall is filled with high grade mineral glass wool.

Temperature range from
ambient to 70degC is controlled by EGO/ JUMO German capillary type thermostat
with an accuracy of +/- 0.5degC.

The unit is provided
with manual tilting device, so as to tilt the tray holding the eggs.

To work on 220/230 volts

A.C. Supply.

Capacity

50 eggs.

100 eggs.

200 eggs.

250 eggs.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003