

Product Name :
Naugra Lab Environmental Stability Chamber

Product Code :
NLAB-LABORATORYEM70003



Description :

Naugra Lab Environmental Stability Chamber

Technical Specification :

Environmental Stability Chamber

Environmental Stability Chamber designed as per ICH Guidelines and to comply requirements of GMP for carrying out various stability and shelf life tests on Pharma, Bulk Drugs, Drug Substances, Biomedical, Healthcare, Food Products, Biotech etc.

The chambers are manufactured as per GMP regulations.

Constructed and fabricated on fully automatic CNC machines.

Environmental Stability Chambers Inner & Outer both made of grade stainless steel.

Outer chamber made of Mild Steel Sheet or grade stainless steel (Option).

PUF insulation on fully automatic High Pressure Puf Injector machine leaving no gaps for higher uniformity and stability.

The unit has full length inner plexi glass door and outside metal door with magnetic gasket & lock.

The unit is supplied with perforated Stainless Steel trays.

Caster wheels for easy moveability.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003