

**Product Name :**  
Naugra Lab Filtration Assembly

**Product Code :**  
NLAB-LABORATORYEM65020



**Description :**

Naugra Lab Filtration Assembly

**Technical Specification :**

Filtration Assembly:-

Complete kit, ready to use. Ideal for filtering HPLC solvents, or any other general purpose filtration.

Kit comprises of:-

All Glass Membrane Filter Holder, Borosilicate Glass, 47mm, with Sintered disc for membrane support and aluminum 'duck' clamp.

Vacuum Pump, Model: P-F6 Compact and efficient.

Silicon Rubber Tubing 6x9mm -1 mtr.

Pack of 100 Membranes, 47mm.

Solvent Filtration Assembly with Pack of Cellulose Nitrate Membranes 0.2 m.

Solvent Filtration Assembly with Pack of Nylon Membranes 0.2 m

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003