

Product Name :
Naugra Lab Flocculator Jar Test Apparatus

Product Code :
NLAB-LABORATORYEM64001



Description :

Naugra Lab Flocculator Jar Test Apparatus

Technical Specification :

Flocculator Jar Test Apparatus

Features:-

Accepts two, four & six samples

Digital indication of speed and time

Two pre-set programs

Built in illumination for easy viewing of samples

Flocculation/Jar Test Apparatus, is a flocculation (jar) tests on water and effluent samples.

A digital display clearly indicates the speed of the rotational stainless steel paddles, which can be varied from 25 to 250rpm.

The digital timer can be set to count down from 1 to 99 minutes.

After count down, the paddles stop and an alert sounds.

The timer then counts up in minutes in order to measure settling time.

Two pre-set programs allow routine speed and time parameters to be input and then selected by simply pressing one button.

Parameters are stored in the memory and can be recalled or altered at any time.

The intuitive touch sensitive keypad is totally flat for easy cleaning.

Specifications:-

No. of Jar	
Speed range	
Pre-set programs	

Time 9 minutes	
Controls Keypad	
Digital Display	
Power, 50Hz, 100W	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003